

GENESY 1500

ELECTROTHERAPY 4 CHANNELS | 414 PROGRAMS

PROFESSIONAL FOUR-CHANNEL PRODUCT USED IN THE MEDICAL FIELD, COMPLETE OF ALL THE CURRENTS FOR PHYSIOTHERAPY AND REHABILITATION THERAPY.

92 ELECTROTHERAPY MENU (EMS/RUSSIAN) | 84 ACTION NOW | 54 SERIAL SEQUENTIAL STIMULATION | 53 SPORT
33 PAIN | 29 FITNESS | 23 MICROCURRENTS | 16 BEAUTY | 11 REHAB | 9 DENERVATED | 6 INTERFERENTIALS
3 INCONTINENCE | 1 IONOPHORESIS

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

4 independent (8 electrodes)

Frequency

0,3-150Hz *

Pulse amplitude

40-450µs *

Power

0-120mA per channel

Power supply

- Mains
- Rechargeable Batteries

Size

mm 100x160x35


Weight

gr 454


* according to the type of current

EQUIPMENT

- 1 Bag
- 1 GENESY 1500 stimulator
- 4 Cables for electrode connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Power supply unit
- 1 Operating manual



CE
0476



TENS



EMS



IONO



3S



ACTION NOW



RUN TIME



LAST 10



FAVORITES



MULTI USER



PROGRAMMABLE



STIM LOCK



FUNCTION 2+2



MY TRAINER



AUTO STIM



SYNCHRO STIM



WORK TIME



MAINS POWER SUPPLY



RECHARGEABLE BATTERIES

GENESY 1500

ELECTROTHERAPY 4 CHANNELS | 414 PROGRAMS

PROGRAMS LIST

REHAB

Swollen ankles

Atrophy recovery

Hemiplegy-upper limbs

Hemiplegy-lower limbs

Recovery after ACL surgery

Functional recovery

Ankle re-education

Leg re-education

Muscle Spasms

Shoulder subluxation prevention

Vastus medialis reinforcement

Motor point pen

Agonist/Antagonist

Muscle reinforcement

PAIN

High frequency antalgic tens

Low frequency antalgic tens

Endorphinic tens

Knee pain

Post-surgery pain

Chronic pain

Shoulder pain

Carpal tunnel

Trapezius pain

Knee osteoarthritis

Rotator cuff tendinitis

Muscle pain

Nerve compression

Spinal compression

Osteoarthritis

Spinal osteoarthritis

Spinal osteoporosis

Hip osteoarthritis

Ankle osteoarthritis

Muscle tendon injury

Post-surgical pain

Sciatica

Chronic lumbago

Cervical pain

Acute pain post inguinal hernia

Bursitis-tendinitis

Muscle injuries

Bone fractures

Whiplash

Menstrual pain

Conventional antalgic tens

Epicondylitis

Modulated antalgic tens

IONOPHORESIS

Ionophoresis

INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

MICROCURRENTS

Epicondylitis

Scapulohumeral periarthritis

Muscle restoration

Contusion

Edema

Skin ulcera

Sciatica

Lumbago

Brachial neuralgia

Acute pain

Articular pain

Stiff neck

Whiplash

Cervical spondylosis

Shoulder sprain

Carpal tunnel

Knee sprain

Knee osteoarthritis

Ankle sprain

Achille tendon inflammation

Patella tendon inflammation

Rotator cuff inflammation

Tendon inflammation

DENERVATED

Triangular 1

Triangular 2

Triangular 3

Trapezoidal 1

Trapezoidal 2

Trapezoidal 3

Rectangular 1

Rectangular 2

Rectangular 3

INTERFERENTIALS

Knee osteoarthritis

Acute lumbago

Cervical pain

Chronic lumbago

Frozen shoulder

Post-surgical pain

SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Aerobic endurance

Reactivity

Post-competition recovery

Decontracting

Hypertrophy

FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

BEAUTY

Drainage

Lipolysis

Toning massage

Skin tone improvement

Connective massage

Post-pregnancy drainage

Post-pregnancy lipolysis

Post-pregnancy firming

SERIAL SEQUENTIAL STIMULATION

The "3S" program list includes 54 parameter combinations.

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2.

The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.

ACTION NOW

The Action Now program list includes 84 parameter combinations.

ELECTROTHERAPY

The Electrotherapy includes 48 combinations of EMS and 44 Russian (KOTZ) parameters.